

Monitor Neurological Status During Any Conscious Procedure

Choose the B-SMART[®]

Buchbinder Sensory Motor Activated Response Timer

- **Monitor neurological function during conscious sedation for any procedure**
- **Monitor level of consciousness during administration of neuroleptic anesthesia**
- **Alerts OR staff when a patient is not responding to a vibratory stimulus**
- **Alerts OR staff to changes in a patient's neurological status simply, reliably, and accurately**
- **Timed repetitive signals can be set for 30, 60 or 120 second intervals**
- **Patient interactive cooperation is an added safety measure**
- **Easy set-up**



GLN: 0894912002005

VASCULAR INSTRUMENTS
KOVEN
TECHNOLOGY INCORPORATED

Koven Technology, Inc.
Your Vascular Healthcare Partner

T: 800.521.8342 / 314.542.2101 • F: 314.542.6020 • info@koven.com • www.koven.com

 Find us on Facebook