

Sensitive to your Doppler Demands

Choose the ES-100VX MiniDoppler® for quick, simple, quality sounds

- **Ideal for quick pedal pulse checks, distal pulses, and fetal heart sounds** (when used with 2.25 MHz fetal probe)
- **Wide selection of interchangeable probes** (available in 2.25 MHz Fetal; 5, 8, and 10 MHz Vascular)
- **Compact & battery operated**
- **Doppler rollcart available for user convenience**



*Convenient, handheld design ...
Quality sounds*



Doppler Rollcart

GLN: 0894912002005

VASCULAR INSTRUMENTS
KOVEN
TECHNOLOGY INCORPORATED

Koven Technology, Inc.
Your Vascular Healthcare Partner



Find us on Facebook

T: 800.521.8342 / 314.542.2101 • F: 314.542.6020 • info@koven.com • www.koven.com