

Sensitive to your Doppler Demands

Choose the *ES-100X Surgical MiniDop*[®]

- Wide selection of interchangeable probe styles and frequencies for user flexibility
- Autoclavable and single-use probes available
- Specialty probes for Vascular, Laparoscopic, Plastic & Reconstructive, Cardiovascular, and Neurosurgery
- Surgical and non-invasive probes available in 5, 8, 10 and 20 MHz
(non-invasive available in 5, 8 & 10 MHz)
- 1-year warranty on Doppler and probes reduces replacement costs
(See warranty information for CRP, NRP, and ACP series probes)
- Compact and battery operated
- Doppler rollcart and footpedal remote available for convenience & surgeon control in the OR



*Compact & Convenient...
Ideal for the OR*



Optional accessories include a Doppler rollcart and footpedal remote

GLN: 0894912002005



Koven Technology, Inc.
Your Vascular Healthcare Partner



Web based training available for unlimited access to on-demand training & support



Pencil Probe
(P8M05S8A 8 MHz)
* also available in 5 & 10 MHz



Single-Use
(VDP-8 8 MHz)



Surgical Specialty Probes
(8, 10, & 20 MHz available)
Cardiovascular, Laparoscopic, Neurovascular, autoclavable & flat monitoring models available