

Sensitive to your Doppler Demands

Choose the HD-307 Bidirectional Surgical Doppler

- **Wide selection of interchangeable probe styles for user flexibility**
- **Autoclavable and single-use probes available**
- **Specialty probes for Vascular, Laparoscopic, Plastic & Reconstructive, Cardiovascular, and Neurosurgery**
- **Surgical probes available in 8, 10 and 20 MHz**
- **Compact and battery operated**
- **1-year warranty on Doppler and probes reduces replacement costs**
(See warranty information for CRP, NRP, and ACP series probes)
- **Single-Frequency Bidirectional Doppler**
(8, 10 & 20 MHz models available)
GTIN: 00894912002104 (8 MHz)
GTIN: 00894912002111 (10 MHz)
GTIN: 00894912002128 (20 MHz)
- **Doppler rollcart and footpedal remote available for convenience & surgeon control in the OR**



Web based training available for unlimited access to on-demand training & support



Doppler Rollcart

Footpedal Remote



Lightweight & Portable ...
Bidirectional Flow Assessment



Surgical Specialty Probes
(8, 10, & 20 MHz available)
Cardiovascular, Laparoscopic, Neurovascular
autoclavable & flat monitoring models
available



Pencil Probe
(PP8M05N5A 8 MHz)
* also available in 10 MHz)



Single-Use
(VDP-8 8 MHz)

GLN: 0894912002005



Koven Technology, Inc.
Your Vascular Healthcare Partner