

Diagnosis is the First Step

Choose the Smartdop® 30EX Vascular Testing System for early diagnosis of vascular disease

Leasing options available




Durable & Effortless...
Automated, reliable
vascular testing

- Fully customizable - choose only the options & accessories to meet the specific needs of your facility
- Built-in automatic cuff inflator - quickly performs automated ABI, PVR and PPG toe pressure studies
- Automatically calculates and prints ABI and TBI results
- Unique, real-time waveform LCD - displays the waveform before printing
- 30 waveform memory storage - save waveforms to view later or download to Smart-V-Link® software for EHR's & DICOM® compatible files
- Automatic gain and baseline control - no manual adjustments required
- Waveform strip generated for immediate documentation for reimbursement
- Optional PPG probe available for expanded testing capabilities
- Flat vascular probe also available - eliminates the need to hold the probe during testing
- Doppler rollcart available for user convenience
- Spanish display option



Web based training available for unlimited access to on-demand training & support



Smart-V-Link® vascular testing software is available for EHR's & DICOM® compatible files without usage fees or significant hardware costs (run on any Windows® based computer).

GLN: 0894912002005

T: 800.521.8342 / 314.542.2101 • F: 314.542.6020 • info@koven.com • www.koven.com  Find us on Facebook

Koven Technology, Inc.
Your Vascular Healthcare Partner