

Sensitive to your Surgical Probe Demands

- **Durable and Cost-effective**
- **Surgical probes are interchangeable with Koven Surgical Dopplers**
- **Specialty probes for Vascular, Cardiovascular, Laparoscopic, Neurosurgery, Plastic and Reconstructive**
- **Pre-sterile single-use & autoclavable probes available**
- **Wide variety of sterilization options available (ETO, Steris®, Steris® V-Pro, Sterrad®, Autoclave)**
- **Non-invasive probes also available (2, 5, 8, and 10 MHz)**
- **1 Year warranty ***

* 1 year parts & labor warranty when used as recommended. Excludes VDP-8 single-use. Warranties for ACP, NRP, and CRP series probes differ. See warranty information for details.

Cardiovascular Probes

Models: CRP-10H & CRP-20H (10 & 20 MHz)
Sterilization: ETO, Sterrad®, Steris®, Steris® V-Pro
*Adaptor required

Bayonet Style Neurovascular Probe

Models: NRP-10H & NRP-20H (10 & 20 MHz)
Also available in **Flexible Tip** models HP-10 & HP-20 (10 & 20 MHz) for easier arterial access to difficult locations
Sterilization: ETO, Sterrad®, Steris®, Steris® V-Pro
*Adaptor required

Laparoscopic Probe

Model: LRP-8 (8 MHz)
Sterilization: ETO, Sterrad®, Steris®, Steris® V-Pro
*Adaptor required

Flat Monitoring Probe

Model: FDP-8 (8 MHz)
Sterilization: ETO, Sterrad®, Steris® V-Pro
*Adaptor required

Air Emboli Flat Monitoring Probe

Model: BF2M20S8A (2 MHz)
Sterilization: Sterrad®
Compatible with ES-100VX and Bidop® 3 Noninvasive Dopplers

Surgical Pencil Probes

Models: P5M05S8A, P8M05S8A, P10M05S8A (5, 8 & 10 MHz)
Models for HD-307: PP8M05N5A, PP10M05N5A (8 & 10 MHz)
Sterilization: ETO, Sterrad®

Autoclavable Probe

Model: ACP-8 (8 MHz)
Sterilization: Autoclave
* Adaptor required

Surgical Thin Probes

Models: VRP-8 & VRP-20 (8 & 20 MHz)
Sterilization: ETO, Sterrad®, Steris®, Steris® V-Pro
*Adaptor required

Single-Use, Pre-Sterile Probes

Model: VDP-8 (8 MHz)
*Adaptor required



Koven Technology, Inc.
Your Vascular Healthcare Partner